

ECOPALZER



ECOPALZER Sterilizer
(Pre Vacuum System)

PS-140R/P



From Japan, to the world
Safety Formaldehyde-gas Sterilizer

PS-140 R/P

ECOPALZER Sterilizer

Revolutionary New Products “Formaldehyde-Gas Sterilizer ECOPALZER” **CHARACTER**

- 1 Short time Processing**
One operation 310 min
- 2 Low Running Cost**
Pharmaceutical expense + Electricity bill = about US \$4
- 3 Harmless to environment**
- 4 Guarantee of Operator Safety**
No-formaldehyde-gas-emission, odorless
- 5 Easy Operation**
By touching screen display according to purpose of use
- 6 Various Subjects under Low Temperature Sterilization**



Medical Device Approved No. 22100 BZX 00829000




D&D system

Decomposition & Detoxification

Sterilized gas emission to the outside will cause environmental pollution. Needless to say; the strong odor and inhalation of residual formaldehyde-gas is harmful to health.

“**D&D system**”, developed with environmental preservation in mind, is a radical system making detoxification and deodorization of formaldehyde-gas possible in a short time.


“**D&D system**” has changed the future of formaldehyde-gas sterilization procedures.




The liquid crystal color touch panel is used for the final controlling element.

The screen display with the interactive system allows for easy operation. Sterilization process, the remaining time and any information are comprehensibly displayed. In addition, the pressure and the temperature in the process are indicated in the graph.

■PS-140R



■PS-140P



■PS-140R/P

Type	Size(widexdepthxheight)mm	Power Supply	Piping
PS-140R	Out 840 x940x1745 In 415x800x415	Three phase 200V /18A	60φ
PS-140P	Out 840 x1012x1745 In 415x800x415		



MEDIATE CO.,LTD

117-5 Megawa,Makishima-cho,Uji Kyoto
611-0041 Japan
Tel:+81-774-28-1170 Fax:+81-774-28-1169

Homepage <http://www.mediate-gp.co.jp>
E-mail ecopal@mediate-gp.co.jp

